

## Anti-IL-6/IFN $\beta$ 2 Antibody

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-IL-6/IFN $\beta$ 2 Antibody is a monoclonal antibody that specifically neutralizes interleukin-6 (IL-6). As a pleiotropic cytokine playing a central role in inflammatory responses and the tumor microenvironment, IL-6's blockade by this antibody effectively inhibits the activation of the JAK/STAT pathway, making it useful for studying rheumatoid arthritis, cytokine storms, and various solid tumors.
Targets(IC50)	Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481