

A2 glycan (G0), procainamide labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	A2 glycan (G0), procainamide labelled (A2 N-linked oligosaccharide, procainamide labelled), is a multifunctional fluorescent linker consisting of a polymeric N-glycoprotein. The resulting conjugate mimics the antenna elements of N-glycans, offering high sensitivity and specificity.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481