

FA2G2 glycan (G2F)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FA2G2 glycan (G2F) [also known as FA2G2 N-linked oligosaccharide; F(6)A2G2] is a complex N-glycoprotein that functions as a versatile fluorescent linker. The resulting conjugates mimic the antenna component of N-glycans, offering high sensitivity and specificity.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481