

A4 glycan, 2-AB labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	A4 glycan, 2-AB labelled (A4 N-linked oligosaccharide, 2-AB labelled) is a polymeric N-glycoprotein serving as a versatile fluorescent linker. The resulting conjugate mimics the antennary structures of N-glycans, offering high sensitivity and specificity.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481