

A3G3 glycan (G3), procainamide labelled

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | A3G3 glycan (G3), procainamide labelled (A3G3 N-linked oligosaccharide, procainamide labelled; A3G(4)3 glycan, procainamide labelled) is a polyglycosylated N-protein that serves as a versatile fluorescent linker. The resulting conjugate exhibits high sensitivity and specificity by mimicking the antennary elements of N-glycans. |
| Targets(IC50) | Others   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481