

GPEP FA2-KVANKT glycopeptide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	GPEP FA2-KVANKT glycopeptide (F(3)A2 glycan-KVANKT & F(6)A2 glycan-KVANKT) is a type of N-glycoprotein that serves as a multifunctional fluorescent linker. The resulting conjugates demonstrate high sensitivity and specificity by simulating the antenna elements of N-glycans.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481