

## A3 glycan, procainamide labelled

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Store at -20°C**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	A3 glycan, procainamide labelled (A3 N-linked oligosaccharide, procainamide labelled), is a polymeric N-glycoprotein functioning as a versatile fluorescent linker. The resulting conjugate exhibits high sensitivity and specificity by mimicking the antennae elements of N-glycans.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481