

Anti-MPL/TPOR/CD110 Antibody (TA136)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-MPL/TPOR/CD110 Antibody (TA136) is a humanized antibody targeting MPL/TPOR/CD110, suitable for use in related research within the life sciences field.
Targets(IC50)	Immunology/Inflammation related
In vitro	In Ba/F3 cells expressing human MPL, Anti-MPL/TPOR/CD110 Antibody (TA136) analogs stimulated cell proliferation and activated the JAK2/STAT signaling pathway [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481