

Anti-IL-3Ra/CD123 Antibody (SNG-CD123A)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-IL-3Ra/CD123 Antibody (SNG-CD123A) is a humanized monoclonal antibody targeting CD123. CD123 is highly expressed in hematological malignancies such as AML. This antibody specifically binds to CD123, blocking its signaling and mediating potent antibody-dependent cellular cytotoxicity (ADCC) to eliminate tumor cells, serving as a critical tool for researching targeted therapies in hematology.
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481