

Anti-Mesothelin Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Anti-Mesothelin Antibody is a monoclonal antibody targeting. Mesothelin is a cell surface glycoprotein anchored by glycosyl phosphatidylinositol (GPI), which is highly expressed in mesothelioma, pancreatic cancer and ovarian cancer, but limited in normal tissues. Anti-Mesothelin Antibody is commonly used in the research of tumor targeted therapy and the development of antibody drug conjugate (ADC) by specifically binding to mesothelin. |
| Targets(IC50) | Mesothelin |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481