

## Anti-SCN11a/Nav1.9 Antibody

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-SCN11a/Nav1.9 Antibody is a humanized monoclonal antibody that binds specifically to the extracellular domain of Nav1.9 (e.g., $\alpha$ subunit) with high affinity and blocks its mediated ion currents, serving as an important pharmacological tool for researching chronic pain, inflammatory pain, and hereditary pain disorders.
Targets(IC50)	Sodium Channel

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481