

Anti-RAMP3 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-RAMP3 Antibody is a highly specialized humanized monoclonal antibody target and neutralize the receptor activity-modifying protein 3 (RAMP3), inhibition of pathological signaling pathways to suppress tumor-induced angiogenesis and metastatic progression across various experimental models .
Targets(IC50)	Others
In vitro	In 4T1 breast cancer cell assays, Anti-RAMP3 Antibody significantly inhibited cell motility and invasive potential in a dose-dependent manner [1].
In vivo	In mouse models of 4T1 breast cancer, systemic administration of Anti-RAMP3 Antibody effectively suppressed primary tumor growth and lung metastasis [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481