

Anti-TNFRSF21/DR6/CD358 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-TNFRSF21/DR6/CD358 Antibody is a humanized monoclonal antibody targeting death receptor 6 (DR6, also known as TNFRSF21). It can be used to study DR6-mediated apoptotic signaling pathways, the mechanisms underlying neurodegenerative diseases, and the regulation of the tumor immune microenvironment.
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481