

Anti-DLL4 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|-------------|--|
| Description | Anti-DLL4 Antibody is a highly specialized monoclonal antibody target the δ -like ligand 4 within the Notch signaling axis, disruption of ligand-receptor engagement to modulate both vascular organization and tumor cell differentiation across various preclinical experimental oncology models. |
| In vivo | In pancreatic cancer xenograft models, Anti-DLL4 Antibody (10-80 mg/kg) inhibited tumor growth and nonproductive angiogenesis while mitigating liver toxicity [1]. |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481