

Anti-CAPRIN1 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-CAPRIN1 Antibody is a human-derived antibody expressed in CHO cells, targeting CAPRIN1. It features a hulgG1 heavy chain and a huk light chain, with an estimated molecular weight (MW) of 150 kDa. For its isotype control, refer to HumanIgG1kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481