

Anti-CXCL4/PF4 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at low temperature

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-CXCL4/PF4 Antibody is a humanized monoclonal antibody that binds specifically to the chemokine CXCL4 (PF4) with high affinity. By blocking its mediated biological effects and receptor interactions, it inhibits fibroblast activation and pro-inflammatory signaling. It serves as an important pharmacological tool for researching systemic sclerosis, thrombotic disorders, and inflammatory immune responses.
Targets(IC50)	CXCR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481