

A2[3]G1 & A2[6]G1 glycan (G1), 2-AA labelled

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	A2[3]G1 & A2[6]G1 glycan (G1), 2-AA labelled, is a type of poly N-glycoprotein and serves as a versatile fluorescent linker. The resulting conjugates, which mimic structural elements of N-glycan antennae, exhibit high sensitivity and specificity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481