

FA2G2S1 glycan (G2FS1), 2-AB labelled

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	The compound FA2G2S1 glycan (G2FS1), which is 2-AB labelled (FA2G2S1 N-linked oligosaccharide, 2-AB labelled), is a type of polymeric N-glycoprotein that serves as a multifunctional fluorescent linker. The resulting conjugate, by mimicking the antenna elements of N-glycans, exhibits high sensitivity and specificity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481