

PDL-241

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PDL-241 is a humanized monoclonal antibody targeting SLAMF7 (also known as CD319 or CS1), which is highly expressed on multiple myeloma cells and natural killer (NK) cells. It exerts anti-tumor effects by directly inducing apoptosis, mediating ADCC effects, and boosting anti-tumor immunity through NK cell activation, making it widely used in research on targeted therapies for hematological malignancies.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481