

FA2[3]G1 & FA2[6]G1 glycan (G1F), APTS labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FA2[3]G1 & FA2[6]G1 glycan (G1F), APTS labelled, is a type of N-glycoprotein with multiple glycans and serves as a multifunctional fluorescent linker. The resulting conjugate, characterized by high sensitivity and specificity, mimics the antenna structures of N-glycans.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481