

A2[3]G1 & A2[6]G1 glycan (G1), 2-AB labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	A2[3]G1 & A2[6]G1 glycan (G1), 2-AB labelled is an N-glycan featuring two asymmetric antennary Lewis X (LeX) structures with a 2-AB (2-aminobenzamide) label. SLeX serves as a ligand for the cell adhesion molecule E-selectin, which is specifically expressed at sites of inflammatory lesions. Designing conjugates with SLeX polysaccharides enables targeted drug delivery to inflamed areas.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481