

## M6 glycan, 2-AB labelled

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	M6 glycan (Man6), 2-AB labelled (Mannose-6 N-linked oligosaccharide, 2-AB labelled; Oligomannose 6 glycan, 2-AB labelled) is a polyvalent N-glycoprotein that serves as a multifunctional fluorescent linker. The resulting conjugate, by mimicking the antenna elements of N-glycans, exhibits high sensitivity and specificity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481