

FA2G2 glycan (G2F), 2-AA labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FA2G2 glycan (G2F), 2-AA labelled (FA2G2 N-linked oligosaccharide, 2-AA labelled; F(6) A2G2, 2-AA labelled) is a type of polymetric N-glycoprotein, serving as a versatile fluorescent linker. The resulting conjugates exhibit high sensitivity and specificity by mimicking the antennal components of N-glycans.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481