

## M3 glycan, 2-AB labelled

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	M3 glycan (Man3), 2-AB labelled (Mannose-3 N-linked oligosaccharide, 2-AB labelled; Oligomannose 3 glycan, 2-AB labelled) is a versatile fluorescent linker that serves as a polymeric N-glycoprotein. The resulting conjugate, with its simulated N-glycan antenna elements, exhibits high sensitivity and specificity.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481