

M5 glycan, procainamide labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	M5 glycan (Man5), procainamide labelled (Mannose-5 N-linked oligosaccharide, procainamide labelled; Oligomannose 5 glycan, procainamide labelled) is a multifunctional fluorescent linker that acts as a poly-N-glycoprotein. The resultant conjugate exhibits high sensitivity and specificity by mimicking the antenna-like elements of N-glycans.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481