

## M9 glycan, procainamide labelled

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	M9 glycan (Man9), procainamide labelled (Mannose-9 N-linked oligosaccharide, procainamide labelled; Oligomannose 9 glycan, procainamide labelled) is a type of poly N-glycoprotein and serves as a versatile fluorescent linker. The resulting conjugates possess high sensitivity and specificity by mimicking the antennae elements of N-glycans.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481