

Anti-IL-20Ra Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-IL-20Ra Antibody is a highly specialized monoclonal antibody target the human interleukin-20 receptor α subunit with high affinity and selectivity, inhibition of ligand-induced signaling through the JAK-STAT pathway across various preclinical experimental models to modulate keratinocyte proliferation and inflammatory gene expression.
Targets(IC50)	Others
In vitro	In reconstituted human epidermis models, Anti-IL-20Ra Antibody was used to characterize receptor subunits involved in IL-20 subfamily cytokine-induced keratinocyte hyperproliferation [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481