

## Anti-CXCR3/GPR9/CD183 Antibody

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Anti-CXCR3/GPR9/CD183 Antibody is a potent monoclonal antibody targeting the chemokine receptor CXCR3, acting by competitively blocking the binding of pro-inflammatory chemokines (such as CXCL10) to the receptor's extracellular domain. The compound is an essential tool for investigating effector T cell recruitment, the pathogenesis of autoimmune diseases, and tumor immune microenvironment remodeling, aiming to inhibit excessive inflammatory cell infiltration by intervening in the CXCR3-mediated signaling axis.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481