

## Anti-CDCP1/CD318 Antibody (38 E11)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Anti-CDCP1/CD318 Antibody (38 E11) is a potent monoclonal antibody targeting the extracellular domain of the transmembrane glycoprotein CDCP1, acting by specifically binding to its CUB1 domain. The compound is an essential tool for investigating Src-dependent signaling, tumor cell survival, and distal metastasis mechanisms, aiming to downregulate downstream oncogenic signaling cascades by blocking CDCP1 phosphorylation or inducing its internalization.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481