

## A2G2 glycan (G2)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	A2G2 glycan (G2) (A2G2 N-linked oligosaccharide; A2G(4)2 glycan) is a type of polymeric N-glycoprotein that serves as a versatile fluorescent linker. The resulting conjugate exhibits high sensitivity and specificity by mimicking the antenna elements of
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481