


TDD

Chemical Properties

CAS No. : 55062-34-5
 Formula: C₁₆H₁₇N₃O₂
 Molecular Weight: 283.33
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	TDD (Tryptophandehydrobutyrine diketopiperazine) is a fungal metabolite that can be isolated from <i>Streptomyces spectabilis</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5295 mL	17.6473 mL	35.2945 mL
5 mM	0.7059 mL	3.5295 mL	7.0589 mL
10 mM	0.3529 mL	1.7647 mL	3.5295 mL
50 mM	0.0706 mL	0.3529 mL	0.7059 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481