

Anti-KID3 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-KID3 Antibody is a specific monoclonal antibody targeting the KID3 antigen, which is identical to CD24. CD24 is overexpressed in various epithelial tumors and is considered a key marker for tumor progression and cancer stem cells. This antibody is widely applied in tumor diagnostics and research focused on developing novel immunotherapies by blocking CD24-mediated cell adhesion and signaling.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481