

## Anti-DLL3 Antibody

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Anti-DLL3 Antibody is a humanized monoclonal antibody targeting Delta-like ligand 3 (DLL3). DLL3 is specifically overexpressed in small cell lung cancer (SCLC) and high-grade neuroendocrine carcinomas, with minimal expression in normal tissues. This antibody binds to the extracellular domain of DLL3 with high affinity and induces rapid internalization. It serves as a critical tool for developing targeted therapies, such as ADCs or bispecific antibodies, aimed at delivering payloads into tumor cells to induce apoptosis and inhibit tumor growth.
Targets(IC50)	Gamma-secretase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481