

PB0822

Chemical Properties

CAS No. :

Formula: C19H19FN4O2

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PB0822, labeled with [18F] fluoride, is suitable for in vivo PET imaging. [18F]PB0822 serves as a PET radioligand. With a cLogP value of 1.54, [18F]PB0822 is a good predictor of blood-brain barrier penetration and can be used for Alzheimer's disease research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481