

155H1

Chemical Properties

CAS No. :

Formula: C79H120FN19O23S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	155H1 (Compound 11) is a stapled peptide that covalently binds to hMcl1 (172-323) with an IC50 of 18 nM.
Targets(IC50)	Bcl-2 Family

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481