

DOTA-PSMA-EB-01

Chemical Properties

CAS No. :

Formula: C87H115N15O28S2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DOTA-PSMA-EB-01 (Compound LNC1003) is a specific inhibitor of PSMA with an IC50 of 10.77 nM. It enhances tumor uptake and retention of 177Lu. DOTA-PSMA-EB-01 is applicable in the synthesis or research of radiolabeled conjugates (RDC).
Targets(IC50)	PSMA

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481