

(Z-IETD)2-Rh 110

Chemical Properties

CAS No. :

Formula: C74H86N10O25

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(Z-IETD)2-Rh 110 is a fluorescent substrate for caspase-8 (λ_{ex} = 488 nm, λ_{em} = 523 nm) and can be used to measure caspase-8 activity and monitor apoptosis.
Targets(IC50)	Caspase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481