

5-FAM-TEGEARGSVILLK(5-TAMRA)K-NH2

Chemical Properties

CAS No. :

Formula: C111H148N22O30

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	5-FAM-TEGEARGSVILLK(5-TAMRA)K-NH2 is a substrate for ADAMTS-4/ADAMTS-5, useful in the detection of ADAMTS-4 or ADAMTS-5 activity.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481