

DHMQ

## Chemical Properties

CAS No. :

Formula: C<sub>29</sub>H<sub>32</sub>ClIN<sub>2</sub>O<sub>2</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DHMQ is a near-infrared fluorescent probe that binds to the propylamine group of monoamine oxidase A (MAO-A). It enables real-time tracking of MAO-A activity within mice and cells through fluorescence imaging.
Targets(IC50)	Monoamine Oxidase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481