

## ALPS1 peptide from ArfGAP1

### Chemical Properties

CAS No. :

Formula: C118H168N30O38S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	ArfGAP1's ALPS1 peptide is a curvature-selective peptide derived from the ArfGAP1 ALPS1 motif. This curvature selectivity can be utilized to develop broad-spectrum antiviral peptides.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481