

Lactate Dehydrogenase 4, Human Liver

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Lactate Dehydrogenase 4 (LDH-4), Human Liver (LDH-4 isoenzyme, 1H3M isoenzyme, ld4 isoenzyme) is a biological material or organic compound suitable for use in life science research.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481