

Thalidomide-Pip-N-boc

Chemical Properties

CAS No. : 2636798-54-2


Formula: C₂₃H₂₈N₄O₆

Molecular Weight: 456.49

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Thalidomide-Pip-N-boc is a conjugate of an E3 ligase ligand and a linker. It is utilized in the synthesis of [PROTAC Degradator LLC0424].
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1906 mL	10.9531 mL	21.9063 mL
5 mM	0.4381 mL	2.1906 mL	4.3813 mL
10 mM	0.2191 mL	1.0953 mL	2.1906 mL
50 mM	0.0438 mL	0.2191 mL	0.4381 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481