


## Loteprednol

## Chemical Properties

CAS No. :	129260-79-3
Formula:	C <sub>21</sub> H <sub>27</sub> ClO <sub>5</sub>
Molecular Weight:	394.89
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Loteprednol is a synthetic corticosteroid classified under glucocorticoids. It possesses anti-inflammatory, anti-allergic, and immunosuppressive properties. Loteprednol is applicable in the study of ocular inflammation.
Targets(IC50)	Glucocorticoid Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5324 mL	12.6618 mL	25.3235 mL
5 mM	0.5065 mL	2.5324 mL	5.0647 mL
10 mM	0.2532 mL	1.2662 mL	2.5324 mL
50 mM	0.0506 mL	0.2532 mL	0.5065 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481