


Pentynoic acid STP ester sodium

Chemical Properties

CAS No. :	1807530-14-8
Formula:	C ₁₁ H ₆ F ₄ NaO ₅ S
Molecular Weight:	348.20
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pentynoic acid STP ester sodium is a click chemistry reagent that features terminal propargyl and STP ester groups. It is suitable for protein modification.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8719 mL	14.3596 mL	28.7191 mL
5 mM	0.5744 mL	2.8719 mL	5.7438 mL
10 mM	0.2872 mL	1.436 mL	2.8719 mL
50 mM	0.0574 mL	0.2872 mL	0.5744 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481