

## EM-12-Alkyne-C6-OMs

## Chemical Properties

CAS No. : 2570254-42-9


Formula: C<sub>22</sub>H<sub>26</sub>N<sub>2</sub>O<sub>6</sub>S

Molecular Weight: 446.52

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	EM-12-Alkyne-C6-OMs is a conjugate of an E3 ligase ligand and linker, consisting of 2-(2,6-Dioxopiperidin-3-yl)phthalimidine (EM-12) and the respective linker. It is utilized in the synthesis of PROTAC.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2395 mL	11.1977 mL	22.3954 mL
5 mM	0.4479 mL	2.2395 mL	4.4791 mL
10 mM	0.224 mL	1.1198 mL	2.2395 mL
50 mM	0.0448 mL	0.224 mL	0.4479 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481