

Highly purified Type XI collagen, human sternal cartilage

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Type XI collagen, highly purified and derived from human sternal cartilage (Human Type XI collagen, immunization grade), serves as an immunization-grade protein capable of inducing specific antibody production against this collagen in animals. Additionally, collagen functions as a hydrolytic substrate for MMP.
Targets(IC50)	MMP

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481