

AADvac 1 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AADvac 1 TFA is an active tau peptide vaccine utilized in Alzheimer's disease research. It is formed by conjugating the regulatory peptide 294KDNKHHVPGGS305, which drives tau oligomerization, with keyhole limpet hemocyanin (KLH), and is formulated with aluminum hydroxide.
Targets(IC50)	Microtubule Associated

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481