

## RAG8 peptide

## Chemical Properties

CAS No. :

Formula: C<sub>56</sub>H<sub>98</sub>N<sub>16</sub>O<sub>11</sub>

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RAG8 peptide is an antagonist of protease-activated receptor 4 (PAR 4), which inhibits late-stage platelet activation and reduces thrombosis without altering hemostasis or increasing bleeding risk. It is applicable in atherosclerosis research.
Targets(IC50)	Protease-activated Receptor

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481