

## Tte UvrD Helicase

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Tte UvrD Helicase is a thermostable UvrD helicase initially isolated from <i>Thermoanaerobacter tengcongensis</i> . This enzyme can unwind DNA duplexes that have 3' or 5' single-stranded DNA tails, as well as DNA substrates with blunt ends.
Targets(IC50)	DNA/RNA Synthesis

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481