

## Lactate Dehydrogenase 3, Human Erythrocytes

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Lactate Dehydrogenase3 (LDH-3), Human Erythrocytes (LDH-3, LD3, LD-3, 2H2M), is a biological material or organic compound suitable for use in life science research.
Targets(IC50)	Dehydrogenase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481